



## Product #6621

### TechNiche Phase Change Cooling **BROW PAD**

Powered by CoolPax™

- Hook and loop attaches easily to hard hat harness (hard hat not included)
- Cool Pax™ Phase Change Cooling inserts maintain a comfortable 58°F (14°C)
- 100% Cotton carrier, Polyurethane Insert
- Includes one (1) brow pad cooler and one (1) 6663 Cool Pax™ Insert
- Size: One
- Color: Blue

POWERED BY

## COOLPAX™

CoolPax™ is our Phase Change Material (PCM) cooling technology. It uses a non-toxic liquid that freezes at 58°F/14°C and remains at that temperature for up to 1 hour.

### HOW TO USE :

1. Place CoolPax™ inserts in freezer or ice water for 35-40 minutes
2. Insert activated CoolPax™ into garment pockets
3. CoolPax™ vest or accessory is ready to wear!



### COOLPAX™ Inserts

Polyurethane Outer Material  
**Phase Change Material:** 58°F  
**Activation:** Ice Water or Refrigerator  
**Activated Duration:** up to 1 hour

**FOLLOW US ON SOCIAL MEDIA**

To Stay Up To Date On The Latest Deals & News!

